NI PXI-2501 Specifications

24-Channel FET Multiplexer/Matrix

このドキュメントには、日本語ページも含まれています。

This document lists specifications for the NI PXI-2501 multiplexer/matrix module. All specifications are subject to change without notice. Visit ni.com/manuals for the most current specifications.

Refer to the NI Switches Help for detailed topology information.



Note The 1-wire 48×1 and the 2-wire 24×1 multiplexer topologies can enable a unity gain amplifier to reduce FET settling time. Refer to the *NI Switches Help* for more information about this solution and for detailed topology and pinout information.



Caution To ensure the specified EMC performance, operate this product only with shielded cables and accessories.



Caution Refer to the *Read Me First: Safety and Electromagnetic Compatibility* document at ni.com/manuals for important safety and compliance information.

About These Specifications

Specifications characterize the warranted performance of the instrument under the stated operating conditions.

Typical Specifications are specifications met by the majority of the instrument under the stated operating conditions. The

performance of the instrument is not warranted. Typical specifications are tested at 23 °C.

All voltages are specified in DC, AC_{pk} , or a combination unless otherwise specified.

Input Characteristics

All input characteristics are DC, AC_{rms}, or a combination, unless otherwise specified.

Maximum switching voltage.....±10 VDC, 7 VAC (channel-to-ground)

Overvoltage protection

Signals CH<023>, COM<03>
Powered on or off±25 VDC
Signals AB<01>
Powered on±25 VDC
Powered off±15 VDC
FET switch on resistance
Typical50 Ω
Maximum at 25 °C85 Ω
Maximum at 85 °C100 Ω
Total path resistance
Total path resistance Channel-to-analog bus
1
Channel-to-analog bus
Channel-to-analog bus Typical1,650 Ω
Channel-to-analog bus $ \begin{array}{ccccccccccccccccccccccccccccccccccc$
Channel-to-analog bus $ \begin{array}{cccc} Typical$
$ \begin{array}{llllllllllllllllllllllllllllllllllll$

Cold-junction sensor.....60 µV

Multiplexer topologies (nonamplified)
Offset voltage......5 μV



Typical bandwidth	
$(50 \Omega \text{ source}, 1 \text{ M}\Omega 25$	pF load)
-3 dB	400 kHz
–10 dB	1 MHz

-10 dB1 MHz
Dynamic Characteristics Maximum scan rate15,000 Channels/s
Settling time (+5 V to -5 V step), measured with a 6-inch AB connector to a PXI multifunction DAQ device 0.012% accuracy
With amplifier8.5 μs Without amplifier9.0 μs
0.006% accuracy With amplifier10 μs Without amplifier11.5 μs
0.0015% accuracy With amplifier16 μs Without amplifier18 μs
Settling time (+5 V to -5 V step), measured with a 3-meter cable to a PXI multifunction DAQ device
0.012% accuracy With amplifier21 μs Without amplifier45 μs
0.006% accuracy With amplifier30 μs
Without amplifier60 μs 0.0015% accuracy



Note Certain applications might require additional time for proper settling. For information about including additional settling time, refer to the NI Switches Help.

Trigger Characteristics

With amplifier80 μs Without amplifier160 µs

Input trigger	
Sources	PXI trigger lines 0–7,
	Front panel
Minimum pulse width	
PXI trigger lines	70 ns
Front panel	500 ns
Output trigger	
Destinations	PXI trigger lines 0–7,
	Front panel
Pulse width	1 μs

Physical Characteristics

Relay type	. FET switch
I/O connector	. 68-pin male SCSI
Power requirement	. 1.5W at 5 V, typical 0.36 W at 12 V, typical, 0.36 W at -12 V, typical
Dimensions $(L \times W \times H)$. 3U, one slot, PXI/cPCI module 21.6 × 2.0 × 13.0 cm (8.5 × 0.8 × 5.1 in.)
Weight	. 1750 g (6.5 oz)
Environment	
Operating temperature	. 0 °C to 50 °C
Storage temperature	. –20 °C to 70 °C
Relative humidity	. 5% to 85% noncondensing
	noncondensing
Pollution Degree	Č
Pollution Degree	.2

Shock and Vibration

Operational shock	. 30 g peak, half-sine,
	11 ms pulse
	(Tested in accordance
	with IEC 60068-2-27.
	Test profile developed
	in accordance with
	MIL-PRF-28800F.)
Random vibration	
Operating	. 5 to 500 Hz, 0.3 g _{rms}

Nonoperating 5 to 500 Hz, 2.4 g_{rms} (Tested in accordance with IEC 60068-2-64. Nonoperating test profile exceeds the requirements of MIL-PRF-28800F, Class 3.)

Accessories

Visit ${\tt ni.com}$ for more information about the following accessories.



Caution You *must* install mating connectors according to local safety codes and standards and according to the specifications provided by the connector manufacturer. You are responsible for verifying safety compliance of third-party connectors and their usage according to the relevant standard(s), including UL and CSA in North America and IEC and VDE in Europe.

Table 1. Accessories Available for the NI PXI-2501

Accessory	Part Number
NI TB-2605 multiplexer terminal block (1-wire 48 × 1 multiplexer) (2-wire 24 × 1 multiplexer) (2-wire dual 12 × 1 multiplexers) (2-wire quad 6 × 1 multiplexers) (4-wire 12 × 1 multiplexer)	777878-01
NI TB-2606 matrix terminal block (2-wire 4×6 matrix)	777879-01
TBX-68S terminal block with cold-junction sensor	777716-01
CB-68LP 68-pin unshielded I/O connector block	777145-01
1 m SH68-68S shielded cable	185262-01
2 m SH68-68S shielded cable	185262-02
5 m SH68-68S shielded cable	185262-05

3




Figure 1. NI PXI-2501 Hardware Diagram

The following figures show the pinouts for various topologies.




Figure 2. NI PXI-2501 1-Wire 48 × 1 Multiplexer Pinout

Figure 3. NI PXI-2501 2-Wire 24 × 1 Multiplexer Pinout




Figure 4. NI PXI-2501 2-Wire Dual 12 × 1 Multiplexer Pinout

Figure 5. NI PXI-2501 2-Wire Quad 6×1 Multiplexer Pinout




Figure 6. NI PXI-2501 4-Wire 12×1 Multiplexer Pinout

Figure 7. NI PXI-2501 2-Wire 4×6 Matrix Pinout

Compliance and Certifications

Safety

This product meets the requirements of the following standards of safety for electrical equipment for measurement, control, and laboratory use:

- IEC 61010-1, EN 61010-1
- UL 61010-1, CSA 61010-1



Note For UL and other safety certifications, refer to the product label or the *Online Product Certification* section.

Electromagnetic Compatibility

This product meets the requirements of the following EMC standards for electrical equipment for measurement, control, and laboratory use:

- EN 61326-1 (IEC 61326-1): Class A emissions; Basic immunity
- EN 55011 (CISPR 11): Group 1, Class A emissions
- AS/NZS CISPR 11: Group 1, Class A emissions
- FCC 47 CFR Part 15B: Class A emissions
- ICES-001: Class A emissions



Note For EMC declarations and certifications, refer to the *Online Product Certification* section.

CE Compliance (E

This product meets the essential requirements of applicable European Directives as follows:

- 2006/95/EC; Low-Voltage Directive (safety)
- 2004/108/EC; Electromagnetic Compatibility Directive (EMC)

Online Product Certification

To obtain product certifications and the Declaration of Conformity (DoC) for this product, visit ni.com/certification, search by model number or product line, and click the appropriate link in the Certification column.

Environmental Management

NI is committed to designing and manufacturing products in an environmentally responsible manner. NI recognizes that eliminating certain hazardous substances from our products is beneficial to the environment and to NI customers.

For additional environmental information, refer to the NI and the Environment Web page at ni.com/environment. This page contains the environmental regulations and directives with which NI complies, as well as other environmental information not included in this document.

Waste Electrical and Electronic Equipment (WEEE)



EU Customers At the end of the product life cycle, all products *must* be sent to a WEEE recycling center. For more information about WEEE recycling centers, National Instruments WEEE initiatives, and compliance with WEEE Directive 2002/96/EC on Waste Electrical and Electronic Equipment, visit ni.com/environment/weee.

电子信息产品污染控制管理办法 (中国 RoHS)



中国客户 National Instruments 符合中国电子信息产品中限制使用某些有害物质指令 (RoHS)。关于National Instruments 中国 RoHS 合规性信息,请登录 ni.com/environment/rohs_china。 (For information about China RoHS compliance, go to ni.com/environment/rohs_china.)

LabVIEW, National Instruments, NI, ni.com, the National Instruments corporate logo, and the Eagle logo are trademarks of National Instruments Corporation. Refer to the *Trademark Information* at ni.com/trademarks for other National Instruments trademarks. Other product and company names mentioned herein are trademarks or trade names of their respective companies. For patents covering National Instruments products/technology, refer to the appropriate location: **Help»Patents** in your software, the patents.txt file on your media, or the *National Instruments Patent Notice* at ni.com/patents.

NI PXI-2501 仕様

24 チャンネル FET マルチプレクサ / マトリクス

このドキュメントには、NI PXI-2501 マルチプレクサ / マトリクスモジュールの仕様が記載されています。すべ ての仕様は事前の通知なしに変更されることがありま す。最新の仕様については、ni.com/manualsを参照 してください。

トポロジ.......単線式 48×1 マルチ プレクサ、2線式 24×1マルチプレク サ、2バンク2線式 12×1マルチプレク サ、4バンク2線式 6×1マルチプレクサ、 4線式 12×1 マルチプ レクサ、2線式4×6

トポロジ情報については、『NIスイッチヘルプ』を参照 してください。



メモ 単線式 48×1 および 2 線式 24×1 マル チプレクサのトポロジは、ユニティゲイン増 幅器を有効にして、FET の整定時間を短縮す ることができます。これらの設定方法、およ びトポロジとピン配列情報については、 『NIスイッチヘルプ』を参照してください。

マトリクス



注意 指定された EMC のパフォーマンスを 確保するには、シールドケーブルおよびアク セサリを必ず使用してください。



注意 安全規格の詳細については、ni.com/ manuals から『はじめにお読みください:安全 対策と電磁両立性について』を参照してくだ さい。

仕様値について

「仕様」は、記載された動作条件下で保証される計測器 の性能を示します。

「標準仕様」は、記載された動作条件下で大多数の計測 器が満たす仕様を示します。計測器の性能は保証されま せん。標準仕様は、23℃で検証が行われています。

すべての電圧は特に注釈のない限り、DC、ACnk もしくはその組み合わせとします。

入力特性

すべての入力特性は特に注釈のない限り、DC、 ACrms、もしくはその組み合わせとします。

最大スイッチ電圧......±10 VDC、7 VAC (チャンネル / グランド間)

圖電圧保護
信号 CH<023>、COM<03>
電源投入および切断時 ±25 VDC
信号 AB<01>
電源投入時 ±25 VDC
電源切断時±15 VDC
ET スイッチオン抵抗
標準50 Ω
最大(25℃時)85Ω
最大 (85 ℃時)100 Ω
『ス抵抗合計
チャンネル / アナログバス間
標準1,650 Ω

最大 1,900 Ω チャンネル / コモン間

冰十	1,700 \$2
最大	2,150 Ω
アンプ付きマルチプレク	サトポロジ
オフセット電圧	1.5 mV
冷接点センサ	60 µV
マルチプレクサトポロジ	(増幅なし)

オフセット電圧5 µV



1 900 0

- Id Mad I Id		
RF 性能特性	出力トリガ	DV4 1 11 18 = 45
標準帯域幅	出力先	PXI トリカライン ○○ 7、フロントパネ
(50 Ω 信号源、1 MΩ、25 pF 負荷)		ル
−3 dB400 kHz	パルス幅	1 µs
-10 dB1 MHz		
動特性	物理特性	
	リレータイプ	FET スイッチ
最大スキャンレート15,000 チャンネル /s	I/O コネクタ	68 ピンオス SCSI
整定時間 (+5 ~ -5 V ステップ)、6 インチ AB コネク	所要電力	15 \
タを使用して PXI マルチファンクション DAQ デバイ スまでを測定	//IX 6/J	0.36 W(12 V 時標
0.012% 確度		準)、0.36 W(-12 V
アンプ付き8.5 µs		時標準)
アンプなし9.0 µs	外形寸法(奥行×幅×高さ)	
0.006% 確度		PXI/cPCI モジュール、 21.6×2.0×13.0 cm
アンプ付き10 µs		$(8.5 \times 0.8 \times 5.1 \text{ in.})$
アンプなし11.5 μs	壬旦	
0.0015% 確度	重量	1/50 g (6.5 OZ)
アンプ付き16 µs	環境	
アンプなし18 μs	動作温度	0 ∼ 50 ℃
整定時間(+5~-5 V ステップ)、3 メートルケーブル を使用して PXI マルチファンクション DAQ デバイス	保管温度	−20 ∼ 70 ℃
までを測定	相対湿度	5 ~ 85%
0.012% 確度		(結露なきこと)
アンプ付き21 µs	汚染度	2
アンプなし45 µs	最大使用高度	2.000 m
0.006% 確度		,,000
アンプ付き30 µs	室内使用のみ。	
アンプなし60 μs	耐衝撃 / 振動	
0.0015% 確度	動作時衝撃	是士 20 <i>a (</i> 坐正改
アンプ付き80 μs	圳作时倒掌	嵌入 30 g(千正弦 波)、11 ms パルス
アンプなし160 μs		(IEC 60068-2-27 に準
メモ アプリケーションによっては、より長		拠して試験済み。
い整定時間が必要な場合があります。追加整		MIL-PRF-28800F に準 拠してテストプロファ
定時間についての情報は、『NI スイッチヘル		拠してデストノロファ イルを確立。)
プ』を参照してください。	ランダム振動	, , , , , , , , , , , , , , , , , , , ,
トリガ特性	ランダム派到 動作時	5 ~ 500 Hz 0 3 a
	非動作時	
入力トリガ	71 245 11 ··· 5	(IEC 60068-2-64 に準
ソースPXI トリガライン 0 ~ 7、フロントパネ		拠して試験済み。
ル		非動作時のテストプロ ファイルは
最小パルス幅		MIL-PRF-28800F
PXI トリガライン70 ns		Class 3 の要件を上回
フロントパネル500 ns		る。)

アクセサリ

以下のアクセサリの詳細については、ni.com/jp を参照してください。



注意 メイトコネクタは、必ず地域の安全コードと基準、およびコネクタの製造元によって提供された規格に従って取り付ける必要があります。他社製コネクタの安全適合指令、また該当する基準(北米では UL および CSA、ヨーロッパでは IEC および VDE を含む)に従った使用方法を確認してください。

表 1 NI PXI-2501 対応の NI アクセサリ

アクセサリ	製品番号
NI TB-2605 マルチプレクサ端子台 (単線式 48×1 マルチプレクサ) (2 線式 24×1 マルチプレクサ) (2 バンク 2 線式 12×1 マルチプレクサ) (4 バンク 2 線式 6×1 マルチプレクサ) (4 線式 12×1 マルチプレクサ)	777878-01
NI TB-2606 マトリクス端子台 (2 線式 4×6 マトリクス)	777879-01
TBX-68S 冷接点センサ付端子台	777716-01
CB-68LP 68 ピン非シールド I/O 端子台	777145-01
SH68-68S シールドケーブル、1 m	185262-01
SH68-68S シールドケーブル、2 m	185262-02
SH68-68S シールドケーブル、5 m	185262-05




図 1 NI PXI-2501 ハードウェアダイアグラム

		1						
CH23	35 1	CH47						
CH22	36 2	CH46						
CH21	37 3	CH45						
CH20	38 4	CH44						
CH19	39 5	CH43						
CH18	40 6	CH42						
EXT_TRIG_IN	41 7	GND						
SCAN_ADV	42 8	予約済み						
COM3+ (CH18-CH23) ¹	43 9	COM3- (CH42-CH48) ¹						
COM2+ (CH12-CH17) ¹	44 10	COM2- (CH36-CH41) ¹						
CH17	45 11	CH41						
CH16	46 12	CH40						
CH15	47 13	CH39						
CH14	48 14	CH38						
CH13	49 15	CH37						
CH12	50 16	CH36						
AB1+1	51 17	AB1-1						
AB0+	52 18	AB0-						
CH11	53 19	CH35						
CH10	54 20	CH34						
CH9	55 21	CH33						
GND	56 22	1_WIRE_LO_REF						
CH8	57 23	CH32						
CH7	58 24	CH31						
CH6	59 25	CH30						
COM1+ (CH6-CH11) ¹	60 26	COM1- (CH30-CH35) ¹						
COM0+	61 27	COM0-						
CH5	62 28	CH29						
CH4	63 29	CH28						
CH3	64 30	CH27						
CH2	65 31	CH26						
CH1	66 32	CH25						
CH0	67 33	CH24						
CJS+1	CJS-1							
1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1								
1単線式モードでは使用しない								

CH23+	35 1	CH23-								
CH22+	36 2	CH22-								
CH21+	37 3	CH21-								
CH20+	38 4	CH20-								
CH19+	39 5	CH19-								
CH18+	40 6	CH18-								
EXT_TRIG_IN	41 7	GND								
SCAN_ADV	42 8	予約済み								
COM3+	43 9	COM3-								
COM2+	44 10	COM2-								
CH17+	45 11	CH17-								
CH16+	46 12	CH16-								
CH15+	47 13	CH15-								
CH14+	48 14	CH14-								
CH13+	49 15	CH13-								
CH12+	50 16	CH12-								
AB1+	51 17	AB1-								
AB0+	52 18	AB0-								
CH11+	53 19	CH11-								
CH10+	54 20	CH10-								
CH9+	55 21	CH9-								
GND	56 22	1_WIRE_LO_REF ¹								
CH8+	57 23	CH8-								
CH7+	58 24	CH7-								
CH6+	59 25	CH6-								
COM1+1	60 26	COM1-1								
COM0+	61 27	COM0-								
CH5+	62 28	CH5-								
CH4+	63 29	CH4-								
CH3+	64 30	CH3-								
CH2+	65 31	CH2-								
CH1+	66 32	CH1-								
CH0+	67 33	CH0-								
CJS0+	68 34	CJS0-								
\sim										
12線式モードでは使用しない										

図2 NI PXI-2501 単線式 48×1 マルチプレクサピン配列 **図3** NI PXI-2501 2 線式 24×1 マルチプレクサピン配列

CH23+	35 1	CH23-	CH23+	35 1	CH23-
CH22+	36 2	CH22-	CH22+	36 2	CH22-
CH21+	37 3	CH21-	CH21+	37 3	CH21-
CH20+	38 4	CH20-	CH20+	38 4	CH20-
CH19+	39 5	CH19-	CH19+	39 5	CH19-
CH18+	40 6	CH18-	CH18+	40 6	CH18-
EXT_TRIG_IN	41 7	GND	EXT_TRIG_IN	41 7	GND
SCAN_ADV	42 8	予約済み	SCAN_ADV	42 8	予約済み
COM3+	43 9	СОМ3-	COM3+	43 9	COM3-
COM2+	44 10	COM2-	COM2+	44 10	COM2-
CH17+	45 11	CH17-	CH17+	45 11	CH17-
CH16+	46 12	CH16-	CH16+	46 12	CH16-
CH15+	47 13	CH15-	CH15+	47 13	CH15-
CH14+	48 14	CH14-	CH14+	48 14	CH14-
CH13+	49 15	CH13-	CH13+	49 15	CH13-
CH12+	50 16	CH12-	CH12+	50 16	CH12-
AB1+	51 17	AB1-	AB1+	51 17	AB1-
AB0+	52 18	AB0-	AB0+	52 18	AB0-
CH11+	53 19	CH11-	CH11+	53 19	CH11-
CH10+	54 20	CH10-	CH10+	54 20	CH10-
CH9+	55 21	CH9-	CH9+	55 21	CH9-
GND	56 22	1_WIRE_LO_REF ¹	GND	56 22	1_WIRE_LO_REF ¹
CH8+	57 23	CH8-	CH8+	57 23	CH8-
CH7+	58 24	CH7-	CH7+	58 24	CH7-
CH6+	59 25	CH6-	CH6+	59 25	CH6-
COM1+1	60 26	COM1-1	COM1+	60 26	COM1-
COM0+	61 27	COM0-	COM0+	61 27	COM0-
CH5+	62 28	CH5-	CH5+	62 28	CH5-
CH4+	63 29	CH4-	CH4+	63 29	CH4-
CH3+	64 30	CH3-	CH3+	64 30	CH3-
CH2+	65 31	CH2-	CH2+	65 31	CH2-
CH1+	66 32	CH1-	CH1+	66 32	CH1-
CH0+	67 33	CH0-	CH0+	67 33	CH0-
CJS0+	68 34	CJS0-	CJS0+	68 34	CJS0-
4					
12線式モードでは使用しない			12線式モードでは使用しない		

図4 NI PXI-2501 2 バンク 2 線式 12×1 マルチプレクサピン配列

図 5 NI PXI-2501 4 バンク 2 線式 6×1 マルチプレクサピン配列

```
35 1
                    CH11B-
                                                       COL5+
                                                              35
                                                                       COL5-
    CH11B+
            36
                                                       COL4+
                                                               36
                                                                  2
                                                                       COL4-
    CH10B+
                2
                    CH10B-
                                                               37
            37
                    CH9B-
                                                       COL3+
                                                                  3
                                                                       COL3-
     CH9B+
                3
                                                               38
                                                                  4
     CH8B+
            38
                4
                    CH8B-
                                                       COL2+
                                                                       COL2-
                                                               39
                                                                  5
     CH7B+
            39
                5
                    CH7B-
                                                       COL1+
                                                                       COL1-
                    CH6B-
                                                       COL0+
                                                               40
                                                                  6
                                                                       COL0-
     CH6B+
            40
                6
                                                 EXT_TRIG_IN
                                                               41
                                                                  7
EXT TRIG IN
            41
                7
                                                                       GND
                    GND
                                                               42
                                                                 8
 SCAN_ADV
            42
                8
                    予約済み
                                                   SCAN_ADV
                                                                       予約済み
                                                               43 9
   COM1B+
            43
                9
                    COM1B-
                                                      ROW3+
                                                                       ROW3-
                                                              44
                                                                  10
                                                                       ROW2-
   COM0B+
            44 10
                    COM0B-
                                                      ROW2+
     CH5B+
            45
                11
                    CH5B-
                                                       COL5+
                                                               45 11
                                                                       COL5-
                                                              46
                                                                  12
                                                                       COL4-
     CH4B+
            46
                12
                    CH4B-
                                                       COL4+
                                                       COL3+
                                                              47
                                                                  13
                                                                       COL3-
     CH3B+
            47
                13
                    CH3B-
     CH2B+
            48
                14
                    CH2B-
                                                       COL2+
                                                               48
                                                                 14
                                                                       COL2-
                                                               49
                                                                  15
     CH1B+
            49 15
                    CH1B-
                                                       COL1+
                                                                       COL1-
                                                       COL0+
                                                              50
                                                                  16
                                                                       COL0-
     CH0B+
            50 16
                    CH0B-
     AB0B+
                    AB0B-
                                                 AB1+ (ROW2+)
                                                               51 17
                                                                       AB1- (ROW2-)
            51
                17
                                                 AB0+ (ROW0+)
                                                                       AB0- (ROW0-)
                                                              52
                                                                  18
     AB0A+
            52
                18
                    AB0A-
    CH11A+
                                                       COL5+
                                                               53
                                                                  19
                                                                       COL5-
            53 | 19
                    CH11A-
                                                              54 20
    CH10A+
            54 20
                    CH10A-
                                                       COL4+
                                                                       COL4-
                                                       COL3+
                                                               55
                                                                  21
                                                                       COL3-
     CH9A+
            55 21
                    CH9A-
                    1_WIRE_LO_REF<sup>1</sup>
                                                        GND
                                                              56
                                                                  22
      GND
            56 22
                                                                       1_WIRE_LO_REF<sup>1</sup>
                                                       COL2+
                                                               57
                                                                  23
                                                                       COL2-
     CH8A+
            57
                23
                    CH8A-
                                                       COL1+
                                                                       COL1-
                                                               58
                                                                  24
     CH7A+
            58
                24
                    CH7A-
                                                       COL0+
                                                               59
                                                                  25
                                                                       COL0-
     CH6A+
                25
                    CH6A-
            59
                                                      ROW1+
                                                               60 26
                                                                       ROW1-
   COM1A+
            60 26
                    COM1A-
                                                      ROW0+
                                                               61 27
                                                                       ROW0-
   COM0A+
            61 27
                    COM0A-
                                                       COL5+
                                                              62
                                                                  28
                                                                       COL5-
     CH5A+
                    CH5A-
            62 28
     CH4A+
                    CH4A-
                                                       COL4+
                                                              63 29
                                                                       COL4-
            63
                29
                                                       COL3+
                                                               64
                                                                  30
                                                                       COL3-
     CH3A+
            64 30
                    СНЗА-
                                                       COL2+
                                                               65
                                                                  31
                                                                       COL2-
     CH2A+
                    CH2A-
            65 31
                                                       COL1+
                                                              66
                                                                  32
                                                                       COL1-
     CH1A+
                    CH1A-
            66
                32
                                                       COL0+
                                                                       COL0-
     CH0A+
            67 33
                    CH0A-
                                                               67 33
                                                       CJS+1
                                                              68 34
                                                                       CJS-1
     CJS+1
            68 34
                    CJS-1
     14線式モードでは使用しない
                                                       1マトリクスモードでは使用しない
```

図 6 NI PXI-2501 4 線式 12×1 マルチプレクサピン配列

図7 NI PXI-2501 2 線式 4×6マトリクスピン配列

認可および準拠

安全性

この製品は、計測、制御、実験に使用される電気装置 に関する以下の規格および安全性の必要条件を満たし ます。

- IEC 61010-1、EN 61010-1
- UL 61010-1、CSA 61010-1



メモ UL およびその他の安全保証については、製品ラベルまたは「オンライン製品認証」セクションを参照してください。

電磁両立性

この製品は、計測、制御、実験に使用される電気装置に関する以下の EMC 規格の必要条件を満たします。

- EN 61326-1 (IEC 61326-1): Class A エミッション、 基本イミュニティ
- EN 55011 (CISPR 11): Group 1、Class A エミッション
- AS/NZS CISPR 11: Group 1、Class A エミッション
- FCC 47 CFR Part 15B: Class A エミッション
- ICES-001: Class A エミッション



メモ EMC 宣言および認証については、 「オンライン製品認証」セクションを参照して ください。

CEマーク準拠((

この製品は、該当する EC 理事会指令による基本的要件に適合しています。

- 2006/95/EC、低電圧指令(安全性)
- 2004/108/EC、電磁両立性指令(EMC)

オンライン製品認証

この製品の製品認証および適合宣言(DOC)を入手するには、ni.com/certificationにアクセスして型番または製品ラインで検索し、保証の欄の該当するリンクをクリックしてください。

環境管理

ナショナルインスツルメンツは、環境に優しい製品の設計および製造に努めています。NIは、製品から特定の有害物質を除外することが、環境および NI のお客様にとって有益であると考えています。

環境の詳細な情報については、ni.com/environment (英語) の NI and the Environment を参照してください。このページには、ナショナルインスツルメンツが準拠する環境規制および指令、およびこのドキュメントに含まれていないその他の環境に関する情報が記載されています。

麻電気電子機器(WEEE)



欧州のお客様へ 製品寿命を過ぎたすべての 製品は、必ず WEEE リサイクルセンターへ送 付してください。WEEE リサイクルセンター およびナショナルインスツルメンツの WEEE への取り組み、および廃電気電子機器の WEEE 指令 2002/96/EC 準拠については、 ni.com/environment/weee (英語)を参照 してください。

电子信息产品污染控制管理办法 (中国 RoHS)



中国客户 National Instruments 符合中国电子信息产品中限制使用某些有害物质指令 (RoHS)。关于 National Instruments 中国 RoHS 合规性信息,请登录 ni.com/environment/rohs_china。 (For information about China RoHS compliance, go to ni.com/environment/rohs_china,)

LabVIEW、National Instruments、NI、ni.com、National Instruments のコーポレートロゴ及びイーグルロゴは、National Instruments Corporation の商標です。その他の National Instruments の商標については、ni.com/trademarks に掲載されている「Trademark Information」をご覧下さい。本文書中に記載されたその他の製品名および企業名は、それぞれの企業の商標または商号です。National Instruments の製品 / 技術を保護する特許については、ソフトウェアで参照できる特許情報(ヘルプー特許情報)、メディアに含まれている patents.xxt ファイル、または「National Instruments Patent Notice」 (ni.com/patents) のうち、該当するリソースから参照してください。